

EBAB03

Portable Breathalyser System and Printer

The EBAB03 and its optional printer meet the very latest European standard EN 15964 for Police breathalyser testing. Meeting this specification ensures that it will remain

accurate and reliable. It has multiple modes of operation, and also has the capability of recalling any of the older tests onto its screen for reviewing or reprinting.



- Multiple Test Modes
- Pass/Fail Testing
- Analysis Test
- Screening of liquid for alcohol

PRINTED EVIDENCE

- The time and date of the test
- The last calibration date
- The subjects result and place, including subject name and signature, operator name and signature, witness name and signature



Richard Webb #WAS John Simons Sal Peter Naude Ra Sannyside

SPECIFICATIONS

Analytical Principal	Fuel cell for fast and accurate analysis of breath - such as acetone and hydrocarbons,
	paint and glue fumes, and non-alcoholic foods and medicines
Measuring Range	0.02 to 2.50 mg/L BrAC
Operating Humidity	10 to 95%, non-condensing
Operating Temperature	-5°C to 40°C
Memory	Up to 3,000 test records
Dimensions	120 x 60 x 30 mm
Weight	Approximately 160 grams (including batteries)



MOUTHPIECE TEST

The result will display the exact The result will be displayed as blood alcohol concentration.



CUP TEST

PASS or FAIL.



SNIFFER MODE

To test a bottle or container for alcohol content.







